
HIRSCHMANN®

rotarus®

The high-class peristaltic pump

HIRSCHMANN®

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X21030100

rotarus®

Continuous dispensing, intelligent control

The rotarus® peristaltic pump series sees Hirschmann taking the leap into the world of continuous dispensing pumps.


An universal planned product family is available with a selection of different motors, varying housing safety classes and intelligent control of delivery volumes, covering a broad spectrum of application areas in the lab or industry.

Hirschmann's liquid handling competence is also evident in the sophisticated functions and innovative details of rotarus®.

50 and 100 Watt motors ensure precise delivery in speed ranges from 0.2 to 500 revolutions per minute. This means that media with a high viscosity can also be accurately dispensed. RFID technology is used to detect the pump head. Configuration data for basic parameters can also be stored in this manner and retrieved at any time. The flow and volume models have automatic blockage detection and tube rupture monitoring.

Overview of details

- Rapid tool-free changing of pump heads
- Automatic pump head and tube detection with RFID technology (rotarus® flow and volume)
- Control panel located directly on pump housing
- External TFT touchscreen with user-friendly menu (standard feature on rotarus® flow and volume)
- Calibratable
- Operation with 24 V DC, outdoor use is also possible with appropriate accessories
- USB and RS232 interface for PC connection

Comfortable use and rapid handling

A control panel for the most important functions is integrated directly into the housing of all rotarus® models. In addition, rotarus® flow and volume are equipped with an external control unit with touchscreen.

The pump head connection technology guarantees rapid handling as a standard feature.

Changing is realised quickly and simply in a single action without using tools. The range encompasses a comprehensive selection of single and multichannel pump heads. Pump heads from other manufacturers can also be fitted without difficulty using an adapter.

intelligent

Fully automatic pump head detection using RFID prevents confusion

interactive

USB and RS232 interface for PC connection

variable

The pump housing can also be attached to a stand or stacked

rapid

Pump heads can be changed with a single action

comfortable

The flow and volume models are comprehensive, but can still be programmed easily via the table module

precise media delivery

Three models, numerous applications

The new rotarus® series is being launched with three versions, each with a choice of two different motors and two unit designs.

An intelligent pump range which always provides the appropriate model for a variety of applications.

rotarus® standard

rotarus® standard is a powerful basic model. The cost-effective, precise delivery pump for small delivery volumes and high viscosities is ideal for flexible applications, thanks to the use of diverse pump heads and tubes. The control panel on the housing enables simple and rapid configuration of the most important parameters, and the speed can be programmed directly on the unit.

rotarus® flow

rotarus® flow is the comfortable delivery pump. With control options for velocity and flow speed, flow is particularly suitable for delivery tasks which require a particular delivery speed calibrated in flow rate (e.g. ml/min). Basic functions can be configured quickly and easily using the control panel on the housing. The external control module with a TFT touch-screen provides clear, self-explanatory menu guidance through further functions. Speed in the case of the flow model can be configured directly on the unit. The flow rate is programmed on the table module.

As with all rotarus® pumps, pump heads on the flow version can also be replaced quickly and simply with a single action. Use of RFID technology enhances user comfort and safety through automatic detection of the pump head and tube. The intelligent functions are supplemented by rupture monitoring and blockage detection on the tube. rotarus® flow also has a calibration function and method memory.

rotarus® volume

rotarus® volume, an extremely comfortable and accurate dispensing pump, has positioned itself right from the outset as the top model in the series. In addition to selection options for speed, flow rate, dispensing duration and total volume, volume also has a range of other intelligent functions. Different pump methods with variable cycles and graduated speeds can also be accurately and reliably controlled. In common with the flow model, the control panel on the housing enables rapid and simple configuration of the basic functions.

The external control module with a TFT touch-screen provides clear, self-explanatory menu guidance through further functions.

intelligent range

Appropriate power for individual use

A choice of two different motors with 50 or 100 watt is available for pump drive. They ensure that the required power is provided to suit the application, flow rate and media.

Series 50

The Series 50 is a powerful precision pump drive for low flow rates or highly viscous media.

Depending on the pump head used, the maximum flow rates.

Series 100

The Series 100 is a rapid pump drive for high flow rates.

Depending on the pump head used, the maximum flow rates.

Housing and interior in perfect combination

Two housing versions with different protection classes are available for rotarus®. The housing with white powder coating conforms to protection class IP54 and is ideally suitable for use in the laboratory. The inox stainless steel design with protection class IP65 meets stricter environmental requirements and is the right choice for use in industrial applications in pharmaceutical, medicinal or food production.

Both housing versions have been developed to conform to the fundamental principles of Hirschmann design. Design and functionality also complement each other in rotarus® to produce units which, thanks to their low space requirement, can be easily integrated in an existing laboratory environment or production system, achieving convincing results during use with their sophisticated technology. The smooth housing surfaces can be cleaned without difficulty.

Numerous versions

Speedier change, more versatile functions

The rotarus® range encompasses a comprehensive selection of single and multichannel pump heads. The pump head connection technology was developed for rapid handling. Changing of heads is realised quickly and simply in a single action without using tools.

Pump heads from other manufacturers can also be fitted without difficulty using an adapter. rotarus® single channel pump heads are easily cascaded (i.e. they can be stacked in an assembly and operated on a single rotarus® pump).

PK 10-16

PK 10-24

HIRSCHMANN

MKF 4-8

MKF 8-4

MKF 12-8

Specially developed accessory range for maximum adaptability

Stacking aid
For stacking 2 or more rotarus®.

Stand holder
The stand clip also facilitates use of rotarus® in a stand assembly.

Control unit holder
For placing the control module on the rotarus®.

Accessories for rotarus®	
Order number	Article designation
953 11 20	rotarus® pedal switch
953 11 21	rotarus® mouse
953 11 61	rotarus® system platform
953 11 62	rotarus® control unit holder
953 11 63	rotarus® stacking aid
953 11 70	rotarus® stand holder

Accessories for pump heads		
Order number	Article designation	Suitable for pump head
953 12 50	rotarus® adapter plate	for MASTERFLEX PK
953 12 51	rotarus® adapter plate	for ISMATEK PK


Pedal switch - so you have both hands free.
An ingenious accessory range ensures that you have complete freedom of action. For example, you can start dispensing with the pedal switch and have both hands free for work on the laboratory table.

System platform
Non-slip depositing surface for rotarus® control unit, inclusive cable guide.

Mouse
For extremely sensitive medium dispensing.

Adapter plate
For use with pump heads from other manufacturers.

Connection cable for the serial interface
Cable for linking units to the serial interface (9-pole) and USB interface of the PC.

Tubes

Tygon® R-3603/E-3603

Tygon® laboratory tube for general use in the lab and industry.

Properties:

highly transparent, temperature resistance from -54 to 74 °C, Shore hardness A 55, autoclavable, excellent chemical resistance to the majority of chemical, inorganic lab products, low particle detachment, due to extremely smooth inner wall, extreme flexibility.

Tygon® LFL

Tygon® Long Flex Life pump tube for general applications in the lab and industry and a long service life.

Properties:

highly transparent, temperature resistance from -54 to 74 °C, Shore hardness A 55, autoclavable, max. bending fatigue strength of all transparent Tygon® tubes, extremely low particle detachment, broad chemical resistance, long service life.

PharMed® BPT

Biocompatible tubes for sensitive pharmaceutical and biotechnological applications.

Properties:

beige-coloured, temperature resistance from -50 to +135 °C, Shore hardness A 64, autoclavable, good general chemical resistance, excellent resistance to acid, alkali and oxidation, impermeable to light and UV radiation.

Tygon® 2001

Ultra-pure tube free of softeners for applications in the food industry and aggressive media.

Properties:

transparent, temperature resistance from -78 to 57 °C, Shore hardness A 69, autoclavable, broad chemical resistance, free of softeners.

Norprene® Chemical

The tube for chemical applications.

Properties:

beige-coloured, temperature resistance from -60 to 74 °C, Shore hardness A 61, autoclavable, very high chemical resistance to acids, bases, salts, ketones, alcohols.

Tygon® 3350

Sanitary silicon tube.

Silicon Peroxide

Silicon tube for industrial use.


Properties:

translucent, temperature resistance from -80 to +200 °C, Shore hardness A 50, platinum-cured silicon tube, ultra smooth inner wall, silicon tube with organic peroxide, non-toxic, non-haemolytic, non-pyrogenic.

Properties:

translucent, temperature resistance from -51 to +238 °C, Shore hardness A 55, autoclavable, silicon tube with organic peroxide, non-toxic, resistant to ozone and sunlight.

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