

Pipette Controllers & Pumps

Levo

Pipette Controller

Features

- Ease of use
- Precise pipetting control
- Robust and lightweight
- Appropriate for most of the plastic and glass pipettes from 0.1 ~100ml
- 3µm replaceable Hydrophobic filter
- Easy cleaning and maintenance



Compress the suction bellows before attaching the pipette

Button to expel residual liquid when using blow-out pipette

Aspirating and dispensing comfortably achieved with a sensitive lever

Levo E

Pipette Pump

Features

- Volume capacity 2mL, 10mL and 25mL
- Colour-coded by volume with green, blue and red
- Thumb wheel guarantees precise operation
- Resistant to acids and alkalis
- Easy for cleaning and maintenance



3µm replaceable Hydrophobic filter

Silicone adapter for 0.1~100ml pipette tip, fully autoclavable



Designed for single-handed operation, Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The feature of side lever enables liquid dispensing to be emptied completely, without a droplet remaining.

