

# Electronic Bottle-Top Dispenser dFlow

## dFlow

### Electronic Bottle-Top Dispenser

#### Features

- Motor operation minimizes the grip force
- Electronic control decreases the repetitive strain injury
- The remote control panel makes the operating process without the manual disturbance
- Two Dispensing modes, dispenser and stepper function options



#### Specifications

Volume range	0.1-99.9mL, Max piston lift is 10mL, resolution is 100 $\mu$ L
Volume accuracy	Dispenser: R = 10mL $\pm$ 30 $\mu$ L; CV = 10ml $\pm$ 10 $\mu$ L Stepper: R = 1ml $\pm$ 6 $\mu$ L; CV = 1ML $\pm$ 9 $\mu$ L
Operating temperature range	10°C~30°C
Quality control	DIN EN ISO 8655
Control type	External control